

**Proven
best-in-class
reliability,
year after year.**



Oxlife INDEPENDENCE[®]

Portable Oxygen Concentrator

Our flagship device, the Oxlife **INDEPENDENCE**, sets the standard for **long operating life, reliability, durability and 24/7 performance.**

Equipped with DNA Technology, it is the **DME choice for hospital discharge** and the hospital-to-home transition, meeting the majority of long-term oxygen therapy requirements.



Continuous flow up to 3 LPM
and 12 pulse settings



Proven 24/7 use, year after year
backed by a five-year warranty




Use as POC & stationary
dual billing (E1390 & E1392)



No tanks to refill
or ongoing tank deliveries



24/7 remote monitoring
with DNA Technology




Long-life, dual battery
operation of up to six hours



Intuitive, easy to use
LCD display and setting controls





DNA TECHNOLOGY®

Remote asset & clinical monitoring

O₂ Concepts' proprietary, device-initiated connectivity solution

GPS-enabled inventory control
and loss prevention technology

24/7 remote monitoring
for device performance

Easily access data from any device
via secure, cloud-based platform

Proactively promote enhanced
patient therapy adherence

Reduce field calls with intuitive
customer support resources

Manage hospital DME closet
inventory with device Geofence



OXLIFE INDEPENDENCE® PRODUCT SPECIFICATIONS¹

Dimensions	Height: 20.2" • Width: 10.9" • Depth: 9.45"
Weight²	Device: 16.7 lbs
Oxygen Flow Rate (Pulse 15-40 BPM)	Continuous: 0.5 - 3 LPM in 0.5 increments Pulse: 0.5 - 6.0 in 0.5 increments (8 - 96 ML)
Battery Duration³	5.75 hours @ 2 pulse at 20 bpm 2.5 hours @ 2 LPM continuous
Battery Recharge Time⁴	2.5 hours
Oxygen Concentration	87% - 95% at sea level at all settings
AC Power	Universal Input: 100 - 240V, 50/60Hz 50 - 60 Hz • 2.0A maximum
DC Power	12 - 15 VDC • 15 amps

¹ O₂ Concepts® reserves the right to revise specifications at any time without notice ² Production weight excludes cart kit

³ Battery duration presented for two-battery system ⁴ Battery charge time increases if device is operating while charging

Smarter, more efficient and more reliable oxygen therapy than ever before

A legacy of design innovations, cutting edge technology and unparalleled support

Relied upon by leading academic and research medical centers, and by hospital groups and homecare providers nationwide

The most cost-effective device & operating model

A simplified, effective hospital/clinic discharge process

The ability to remotely monitor and locate every device in your fleet

Designed for the reliability and durability needs of the DME industry

Contact us at **(877) 867-4008 (prompt 6)** or at **info@o2-concepts.com** for more information or to consult with one of our territory managers.

o2-concepts.com



O₂ Concepts®, Oxlife **INDEPENDENCE**®, DNA Technology® and Energy Smart Technology® are trademarks of O₂ Concepts®, LLC and may not be reproduced without express consent.

© 2024 O₂ Concepts®, LLC. All rights reserved.

900-1019 Rev A