



## A leader in the tankless oxygen therapy revolution, the Oxlife **INDEPENDENCE**® is unlike any other POC.

Reliable, durable and easily maneuverable, the Oxlife **INDEPENDENCE**® delivers the continuous flow oxygen therapy you need so you can continue to live the life you want.

Powered by our patented Energy Smart Technology,® the Oxlife **INDEPENDENCE**® provides more than 3x the oxygen flow of other leading portable oxygen concentrators (POCs), the longest battery life available and proven, continuous reliable performance for 24/7 therapy — year after year.



If your healthcare provider determines oxygen therapy is right for you, contact your local provider to ask for the Oxlife **INDEPENDENCE**® by O2 Concepts®.

Oxygen therapy is typically covered by most insurance companies. Check with your insurance company regarding coverage and payment details.

Learn more about the Oxlife **INDEPENDENCE**® by O2 Concepts® at [O2Concepts.com](https://O2Concepts.com).



All O2 Concepts® POCs are proudly innovated, engineered and assembled in the United States.



## Effective oxygen therapy anytime, anywhere — no tank required.


### Oxlife **INDEPENDENCE**®



Meets the majority of long-term oxygen therapy patient requirements



Continuous flow up to 3 LPM



Long-life, dual battery operation of up to six hours



Five-year warranty



Proven 24/7 use as portable and stationary



No tanks to refill or the limitations of ongoing tank deliveries



## Key Features

Continuous flow (up to 3.0 LPM) and pulse flow (.5 to 6) settings

Designed for mobility with an easily maneuverable integrated cart and handle

Long-life, dual battery configuration supports battery swap and charging while device is in use

AC and DC operation in all settings and meets FAA requirements for air travel

Durable aluminum/polycarbonate case with best-in-class components

Compatible for use in conjunction with CPAP and BiPAP



ENERGY SMART TECHNOLOGY®

## Mobility at home, at work and on the go

O<sub>2</sub> Concepts® POCs are conveniently portable around the house, at the workplace and on the road, providing you the mobility you need to maintain your lifestyle.

## No tank headaches required

No ongoing deliveries to coordinate. No bulky tanks to manage. No portable tanks to refill. Just reliable, continuous 24/7 oxygen you can bring along wherever life takes you.

## Stay connected with support

Proprietary DNA® Technology provides performance monitoring and ongoing support to help avoid exacerbations, complications and costly readmissions.